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
Document Title	Part Number
Protocols and User Guides	
C1 System Site Requirements Guide	100-5201
C1 Script Builder Software User Guide	100-7913
Using the C1 System to Capture Cells from Cell Culture and Perform Pre-amplification Using Delta Gene™ Assays Protocol	100-4904
Using the C1 System to Capture Cells from Cell Culture and Perform Pre-amplification Using TaqMan® Assays Protocol	100-6117
Using the C1 System to Capture Cells from Cell Culture and Perform Pre-amplification of MicroRNA Species with TaqMan Assays Protocol	100-6667
Using the C1 System to Generate Single-Cell Libraries for DNA Sequencing Protocol	100-7135
Using the C1 System to Generate Single-Cell cDNA Libraries for mRNA Sequencing Protocol	100-7168
Generate cDNA Libraries with the C1 Single-Cell mRNA Seq HT IFC and Reagent Kit v2 Protocol	101-4964
Fluidigm Single-Cell Preparation Guide	100-7697
Quick References, Templates, and Specification Sheets	
C1 System Specification Sheet	100-5476
C1 System Usage Quick Reference	100-5003
Updating the C1 System Software Quick Reference	100-6217
C1 mRNA Sequencing Data Sheet	101-3386
Single-Cell mRNA Seq PicoGreen Template	100-6260
Using a Microscope with an Automated Stage Quick Reference	100-6130
Minimum Specifications for Imaging Cells in Fluidigm Integrated Fluidic Circuits Specification Sheet	100-5004
Pipetting Maps for the C1 mRNA Seq HT v2 Protocol: Quick Reference	101-4965

[For technical support visit fluidigm.com/support](https://www.fluidigm.com/support)